



March 2017

PELLE CURTENAL MEL

CHARACTERISTIC

Melamine Based Resin

TECHNICAL FEATURES

Chemical Structure : Modified melamine formaldehyde resin
Appearance : Cream colored powder
Concentration (%) : Min. 95,0
pH (10 % in water) : 7,0 - 9,0

PROPERTIES AND APPLICATION

PELLE CURTENAL MEL is a product of melamine formaldehyde resin which is used to fill loosely structured parts of the skin. It improves fullness. It gives fine and smooth grain. It can be used alone or in combination with anionic retanning agents. The best results are obtained by using vegetable extracts and synthetic tannins. It should be added to drum before retanning agents and should be turned for 30 – 45 minutes in any kind of chrome leather.

It is used amount of between 1,0 % and 6,0 % according to the desired leather type.

STORAGE

In a cool and dry place in original packaging, can be stored for a period of one year. Product must be protected from excessive heat and cold.

PACKAGING

25 kg bag.

“TOUCH WITH PELLECURTEX”

Technical information and recommendations contained herein, is reliable, but conditions and methods of use of our products are beyond our control not imply any warranty. We recommend to use to our customer after identified our products and the suitability of proposals of use. Tips on the use of our products be interpreted in terms of encouraging infringement of patent rights of others.